

# Lovibond® Colour Measurement Tintometer® Group



## Pfund

Part Number: 139372



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Pfund

A single filter for checking the PFXi instruments in Pfund scale.

### Technical Data

Glass Standard Values	Approximate value
	- - 30

Filters	1
---------	---

### Delivery Scope

- Pfund 30mm filter
- Certificate of conformity