

Lovibond[®] Colour Measurement Tintometer[®] Group



Lovibond Red, Yellow, Neutral Filter Set (Model F, E, D)

Part Number: 185000



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Lovibond Red, Yellow, Neutral Filter Set (Model F, E, D)

A set of five glass conformance filters for checking a Model F, E and D in Lovibond RYBN scale.

Technical Data

Glass Standard Values	Approximate Values
	- 0.1R 0.5Y
	- 1.0R 2.6Y
	- 1.7R 5.5Y
	- 3.4R 15.0Y
	- 5.0R 34.0Y

Filters

5

Delivery Scope

- Lovibond RYBN 0.1R 0.5Y filter
- Lovibond RYBN 1.0R 2.6Y filter
- Lovibond RYBN 1.7R 5.5Y filter
- Lovibond RYBN 3.4R 15.0Y filter
- Lovibond RYBN 5.0R 34.0Y filter
- Certificate of conformity