

Lovibond[®] Colour Measurement Tintometer[®] Group



Lovibond Red, Yellow, Blue Filter Set (Model F - BS684)

Part Number: 187000



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Lovibond Red, Yellow, Blue Filter Set (Model F - BS684)

A set of ten glass conformance filters for checking a Model F (BS684) in Lovibond RYBN scale.

Technical Data

Glass Standard Values	Approximate Values
	- 0.6R 1.0Y 0.6B
	- 1.1R 2.1Y 0.8B
	- 1.4R 3.4Y 0.8B
	- 1.9R 4.3Y 1.0B
	- 2.3R 7.3Y 1.0B
	- 2.9R 9.7Y 1.1B
	- 4.3R 17.1Y 1.0B
	- 4.7R 19.0Y 0.9B
	- 5.6R 31.0Y 1.0B
	- 7.3R 51.0Y 1.0B

Delivery Scope

- Lovibond RYBN 0.6R 1.0Y 0.6B filter
- Lovibond RYBN 1.1R 2.1Y 0.8B filter
- Lovibond RYBN 1.4R 3.4Y 0.8B filter
- Lovibond RYBN 1.9R 4.3Y 1.0B filter
- Lovibond RYBN 2.3R 7.3Y 1.0B filter
- Lovibond RYBN 2.9R 9.7Y 1.1B filter
- Lovibond RYBN 4.3R 17.1Y 1.0B filter
- Lovibond RYBN 4.7R 19.0Y 0.9B filter
- Lovibond RYBN 5.6R 31.0Y 1.0B filter
- Lovibond RYBN 7.3R 51.0Y 1.0B filter
- Certificate of conformity

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Germany

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
UK

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebu Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaysia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brazil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
USA

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spain

Technical changes without notice

Printed in Germany

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies