

Lovibond® Water Testing Tintometer® Group



DI 10 Dipslide Incubator

For controlled bacterial growth



In order to compare the results of Dipslides and other microbiological tests such as the Coliform / E-coli Test Kit from week to week, the incubation time and temperature must be the same each time the test is performed. This ensures that bacterial growth is always controlled and consistent. The DI10 incubator is therefore an important and necessary step in the analysis of microbiological activity in water. The Lovibond® DI 10 incubator provides reliable incubation in the field, in the laboratory or even on the road. For incubation of dipslides we recommend the purchase of the dipslide holder. This provides a secure support for 12 dipslides in the DI 10 during the incubation period.

- Robust design
- Holds up to 12 dipslides or 10 quanti-discs
- Excellent temperature stability
- In-car operation

Part Number: 56B000701

Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Marine Industry | Oil and Gas | Pharmaceutical Industry | Power and Energy

Technical Data

Display	1 x 16 bit Backlit LCD
Power Supply	12 V 110 V 240 V
Dimensions	246 x 215 x 162 mm

Delivery Scope

- DI 10
- Power supply
- Power cable set (UK, EU and USA)
- Cigarette lighter cable & 8 plug-in adapters
- User manual

Accessories

Title	Part Number
Power cable set	136300
Power supply, 12 V / 48 W	19803550
Dipslide Holder	56B000801

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Germany

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
UK

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaysia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brazil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
USA

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spain